

## REFERENCES

- Adamkin D, Carlo W, Dreyer G, et al: Thermoregulation. In *Neonatal clinical management guidelines*, ed 5, Upper Saddle River, NJ, 2007, Paradigm Health.
- Altmier L: Thermoregulation: what's new? What's not? *Newborn Infant Nurs Rev* 12: 51, 2012.
- Anderson GC, Moore E, Hepworth J, et al: Early skin-to-skin contact for mothers and their healthy newborn infants, *Cochrane Database Syst Rev* 2:CD003519, 2003.
- American Academy of Pediatrics: "Late-preterm" infants: a population at risk, *Pediatrics* 120:1390, 2007.
- American Academy of Pediatrics and American College of Obstetricians and Gynecologists: *Guidelines for perinatal care*, ed 7, Elk Grove Village, Ill, 2012, The Academy.
- Baumgart S: Incubation of the human newborn infant. In Pomerance J, Richardson CJ, editors: *Issues in clinical neonatology*, Norwalk, Conn, 1992, Appleton & Lange.
- Berger I, Marom R, Mimouni F, et al: Weight at weaning of preterm infants from incubator to bassinets: a randomized controlled trial, *Am J Perinatol* 31:535, 2014.
- Bergman N, Linley LL, Fawcus SR: Randomized controlled trial of skin-to-skin contact from birth versus conventional incubator for physiological stabilization in 1200–2199 gram newborns, *Acta Paediatr* 93:779, 2004.
- Bergstrom A, Okong P, Ransjo-Arvidson AB: Immediate maternal thermal response to skin-to-skin care of newborn, *Acta Paediatr* 96:655, 2007.
- Bhatt DR, White R, Martin G, et al: Transitional hypothermia in preterm newborns, *Adv Neonatal Care* 10:515, 2010.
- Billimoria Z, Chawla S, Bajaj M, et al: Improving admission temperature in extremely low birth weight infants: a hospital-based multi-intervention quality improvement project, *J Perinat Med* 41:455, 2013.
- Bissinger RL, Annibale DJ: Thermoregulation in very-low-birth-weight infants during the golden hour, *Adv Neonatal Care* 10:230, 2010.
- Blackburn S, DePaul D, Loan LA, et al: Neonatal thermal care, part III: the effect of infant position and temperature probe placement, *Neonatal Newt* 20:25, 2001.
- Bryanton J, Walsh D, Barrett M, et al: Tub bathing versus traditional sponge bathing for the newborn, *J Obstet Gynecol Neonat Nurs* 33(6):704–712, 2004.
- Cardona-Torres LM, Amador-Licona N, Garcia-Campos ML, Guizar-Mendoza JM: Polyethylene wrap for thermoregulation in the preterm infant: a randomized trial, *Indian Pediatr* 49:129, 2012.
- Castrodale V, Rinehart S: The golden hour: improving the stabilization of the very low birth-weight infant, *Adv Neonatal Care* 14:9, 2014.
- Charpak N, Ruiz JG, Zupan J, et al: Kangaroo mother-care: 25 years after, *Acta Paediatr* 94:514, 2005.
- Cone TE: *History of the care and feedings of the premature infant*, Boston, 1985, Little, Brown.
- Cramer K, Wiebe N, Hartling L, et al: Heat loss prevention: a systematic review of occlusive skin wrap for premature neonates, *J Perinatol* 25:763, 2005.
- DeGuines C, Degrugilliers L, Ghyselen L, et al: Impact of nursing care on the temperature environment in preterm newborns in closed convective incubators, *Acta Paediatr* 102:e96, 2013.
- Dollberg S, Rimon A, Atherton HD, Hoath SB: Continuous measurement of core body temperature in preterm infants, *Am J Perinatol* 17:257, 2000.
- Dollberg S, Mimouni FB, Weintraub V: Energy expenditure in infants weaned from a convective incubator, *Am J Perinatol* 21: 253–256, 2004.
- Duran R, Vatansever U, Acunas B, Sut N: Comparison of temporal artery, mid-forehead skin and axillary temperature recordings in preterm infants <1500 g of birthweight, *J Paediatr Child Health* 45:444, 2009.
- Flenady VJ, Woodgate PG: Radiant warmers versus incubators for regulating body temperature in newborn infants, *Cochrane Database Syst Rev* 4:CD000435, 2003.
- Friedrichs J, Staffileno BA, Fogg L, et al: Axillary temperatures in full-term newborn infants, *Adv Neonatal Care* 13:361, 2013.
- Gaylord MS, Wright K, Lorch K, et al: Improved fluid management utilizing humidified incubators in extremely low birth weight infants, *J Perinatol* 21:438, 2001.
- Godfrey K, Nativio D, Bender CV, Schlenk EA: Occlusive bags to prevent hypothermia in premature infants: a quality improvement initiative, *Adv Neonatal Care* 13:311, 2013.
- Gouchon S, Gregori D, Picotto A, et al: Skin-to-skin contact after cesarean delivery: an experimental study, *Nurs Res* 59:78, 2010.
- Greenspan JS, Cullen AB, Touch SM, et al: Thermal stability and transition studies with a hybrid warming device in neonates, *J Perinatol* 21:167, 2001.
- Haddad L, Smith S, Phillips KD, Heidel RE: Comparison of temporal artery and axillary temperatures in healthy newborns, *J Obstet Gynecol Neonatal Nurs* 41:383, 2012.
- Ibrahim CP, Yoxall CW: Use of plastic bags to prevent hypothermia at birth in premature infants: do they work at lower gestations? *Acta Paediatr* 98:256, 2009.
- Ibrahim CP, Yoxall CW: Use of self-heating gel mattresses eliminates admission hypothermia in infants born below 28 weeks' gestation, *Eur J Pediatr* 169:795, 2010.
- Jia YS, Lin ZL, Lv H, et al: Effect of delivery room temperature on the admission temperature of premature infants: a randomized controlled trial, *J Perinatol* 33:264, 2013.
- Karlsson V, Heinemann AB, Sjors G, et al: Early skin-to-skin care in extremely preterm infants: thermal balance and care environment, *J Pediatr* 161:422, 2012.
- Kim SM, Lee EY, Chen J, Ringer SA: Improved care and growth outcomes by using hybrid humidified incubators in very preterm infants, *Pediatrics* 125:e137, 2010.
- Knobel RB, Wimmer JE, Holbert D: Heat loss prevention for preterm infants in the delivery room, *J Perinatol* 25:304, 2005.
- Knobel RB, Vohra S, Lehman CU: Heat loss prevention in the delivery room for preterm infants: a national survey of newborn intensive care units, *J Perinatol* 25:514, 2005.
- Knobel R, Holditch-Davis D: Thermoregulation and heat loss prevention after birth and during neonatal intensive care unit stabilization of extremely low birth weight infants, *J Obstet Gynecol Neonatal Nurs* 36:280, 2007.
- Leadford AE, Warren JB, Manasyan A, et al: Plastic bags for the prevention of hypothermia in preterm and low birth weight infants, *Pediatrics* 132:e128, 2013.
- Lee GL, Flannery-Bergey D, Randall-Rollins K, et al: Accuracy of temporal artery thermometry in neonatal intensive care infants, *Adv Neonatal Care* 11:62, 2011.
- L'Herault J, Petroff L, Jaffrey J: The effectiveness of a thermal mattress in stabilizing and maintaining body temperature during the transport of very low birth weight newborns, *Appl Nurs Res* 14:210, 2001.
- Lenclen R, Mazarani M, Jugie M, et al: Reducing heat loss in preterm infants at delivery with polyethylene bags, *Arch Pediatr* 9:238, 2002.

43. Lewis DA, Sanders LP, Brockopp DY: The effect of three nursing interventions on thermoregulation in low birth weight infants, *Neonatal Netw* 30:160, 2011.
44. Ludington-Hoe SM, Lewis T, Morgan K, et al: Breast and infant temperatures with twins during shared kangaroo care, *J Obstet Gynecol Neonatal Nurs* 35:223, 2006.
45. Ludington-Hoe SM, Morgan K, Abouelfetoh A: A clinical guideline for implementation of kangaroo care with premature infants of 30 or more weeks' postmenstrual age, *Adv Neonatal Care* 8:S3, 2008.
46. Maayan-Metzger A, Yosipovitch G, Hadad E, et al: Effect of radiant warmer on transepidermal water loss (TEWL) and skin hydration in preterm infants, *J Perinatol* 24:372, 2004.
47. Mathew B, Lakshminrusimha S, Cominsky K, et al: Vinyl bags prevent hypothermia in preterm infants, *Indian J Pediatr* 74:249, 2007.
48. McCall EM, Alderdice F, Halliday HL, et al: Interventions to prevent hypothermia at birth in preterm and/or low birthweight infants, *Cochrane Database Syst Rev* 3:CD004210, 2010.
49. McCarthy LK, Molloy EJ, Twomey AR, et al: A randomized trial of exothermic mattresses for preterm newborns in polyethylene bags, *Pediatrics* 132:e135, 2013.
50. Mellien A: Incubators versus mothers' arms: body temperature conservation in very-low-birth-weight premature infants, *J Obstet Gynecol Neonatal Nurs* 30:157, 2001.
51. Meyer MP: Swaddling and heat loss, *Arch Dis Child Fetal Neonatal Ed* 88:F256, 2003.
52. Mori R, Khanna R, Pledge D, Nakayama T: Meta-analysis of physiological effects of skin-to-skin contact for newborns and mothers, *Pediatr Int* 52:161, 2010.
53. Motil KJ, Blackburn MG, Pleasure JR, et al: The effects of four different radiant warmer temperature set-points used for re-warming neonates, *J Pediatr* 84:546, 1974.
54. Nanndran V, Hoath SB: Thermal management of the low birth weight infant: a cornerstone of neonatology, *J Pediatr* 134:529, 1999.
55. New K, Flenady V, Davies MW: Transfer of preterm infants from incubators to open cot at lower versus higher body weight, *Cochrane Database Syst Rev* 9:CD004214, 2011.
56. Newton T, Watkinson M: Preventing hypothermia at birth in preterm babies: at a cost of overheating some? *Arch Dis Child Fetal Neonatal Ed* 88:F256, 2003.
57. Park HK, Choi BS, Lee SJ, et al: Practical application of kangaroo mother care in preterm infants: clinical characteristics and safety of kangaroo mother care, *J Perinatal Med* 42:239, 2014.
58. Philip AG: The evolution of neonatology, *Pediatr Res* 58:799, 2005.
59. Pinheiro JMB, Boynton S, Furdon SA, et al: Use of chemical warming packs during resuscitation is associated with decreased rates of hypothermia in very low-birth-weight neonates, *Adv Neonatal Care* 11:357, 2011.
- 59a. Russo A, McCready M, Torres L, et al: Reducing hypothermia in preterm infants following delivery, *Pediatrics* 133:e1055, 2014.
60. Sarman I, Can G, Tunell R: Providing warmth for preterm babies by a heated water filled mattress, *Arch Dis Child* 64:29, 1989.
61. Schafer D, Boogaart S, Johnson L, et al: Comparison of neonatal skin sensor temperatures with axillary temperature, *Adv Neonatal Care* 14:52, 2014.
62. Sheldon B: An encapsulated history of thermoregulation in the neonate, *Neo Rev* 5:78, 2004.
63. Silverman W, Sinclair J, Agate F: The oxygen cost of minor changes in heat balance of small newborn infants, *Acta Paediatr Scand* 55:294, 1966.
64. Silverman W, Agate F: Variation in cold resistance among small newborn animals, *Biol Neonate* 6:113, 1964.
65. Simon P, Dunaway D, Bright B, et al: Thermal defense of extremely low gestational age newborns during resuscitation: exothermic mattresses vs polyethylene wrap, *J Perinatol* 31:33, 2011.
66. Sinclair JC: Servo-control for maintaining abdominal skin temperature at 36° C in low birth weight infants, *Cochrane Database Syst Rev* 1:CD001074, 2002.
67. Singh A, Duckett J, Newton T, Watkinson M: Improving neonatal unit admission temperatures in preterm babies: exothermic mattresses, polyethylene bags or traditional approach? *J Perinatol* 30:45, 2010.
68. Smith J, Alcock G, Usher K: Temperature measurement in the preterm and term neonate: a review of the literature, *Neonatal Netw* 32:16, 2013.
69. Smith J, Usher K, Alcock G, Buettner P: Application of plastic wrap to improve temperatures in infants born less than 30 weeks' gestation: a randomized controlled trial, *Neonatal Netw* 32:235, 2013.
70. Thomas KA: Comparability of infant abdominal skin and axillary temperatures, *Newborn Infant Nurs Rev* 3:173, 2003.
71. Thomas KA, Magbalot A, Shinabarger K, et al: Seasonal mapping of NICU temperature, *Adv Neonatal Care* 10:83, 2010.
72. Toubas PL, Nelson R: The role of the French midwives in establishing the first special care units for sick newborns, *J Perinatal* 22:75, 2002.
73. Trevisanuto D, Coretti I, Doglioni N, et al: Effective temperature under radiant warmer: does the device make a difference? *Resuscitation* 82:720, 2011.
74. Trevisanuto D, Doglioni N, Cavallin F, et al: Heat loss prevention in very preterm infants in delivery rooms: a prospective, randomized, controlled trial of polyethylene caps, *J Pediatr* 156:914, 2010.
75. Van der Spek RDG, van Lingen RA, van Zoeren-Grobden D: Body temperature measurement in VLBW infants by continuous skin measurement is a good or even better alternative than continuous rectal measurement, *Acta Paediatr* 98:282, 2009.
76. Vohra S, Roberts R, Zhang B, et al: Heat loss prevention (HELP) in the delivery room: a randomized controlled trial of polyethylene occlusive skin wrapping in very preterm infants, *J Pediatr* 145:750, 2004.
77. Watkinson M: Temperature control of premature infants in the delivery room, *Clin Perinatol* 33:43, 2006.
78. World Health Organization: Thermal protection of the newborn: a practical guide, 1997: [www.who.int/reproductive-health/publications/MSM\\_97\\_2\\_Thermal\\_protection\\_of\\_the\\_newborn/index.htm](http://www.who.int/reproductive-health/publications/MSM_97_2_Thermal_protection_of_the_newborn/index.htm). Accessed March 16, 2014.
79. Zecca E, Corsello M, Priolo F, et al: Early weaning from incubator and early discharge of preterm infants: randomized clinical trial, *Pediatrics* 126:e651, 2010.

---

## RESOURCES FOR PROFESSIONALS

Bissinger RL: Teaching toolbox: method for wrapping the infant in occlusive wrap at birth, *Adv Neonatal Care* 10:239, 2010.